

CCT

COLD COMPRESSION THERAPY

The CCT Professional range of Cold Compression Supports designed to provide rapid relief from injury



Cold Comfort.

Talarmade's effective Cold Compression Therapy system incorporates a gel cold pack covered with frost free lining to allow application directly to the skin. The unique adjustable compression chamber design positioned behind the cold pack provides both contouring and compression around the site of the injury.


Ideal for sports injuries, sprains and strains, ligament damage, general overuse injuries. The Cold Compression Therapy product range comes in eight designs which ensure all body segments can be treated effectively.

Suitable for the treatment of:

- Sports injuries
- Sprains and strains
- Ligament damage
- Overuse injuries

One size fits most.

8 DESIGNS
to ensure
all body segments
can be treated
effectively



The combination of cold and compression drives cold deep into the injured tissue, speeding up the rehabilitation process and aiding faster recovery.

- CCT delivers both cold therapy and compression therapy simultaneously!
- CCT provides an effective way to control the pain and swelling caused by sprains, pulled muscles, injured ligaments and strained tendons.

For best results refrigerate for 2 hours before use and apply to injured area for 10 - 15 minute periods.

Sports Injury Professionals advise
R.I.C.E

REST
ICE
COMPRESSION
ELEVATION



CCT

COLD COMPRESSION THERAPY

CCT provides cold compression support for each area vulnerable to injury. For sports club use convenience products can be purchased as full set which includes one of each CCT support. When purchased as a full set you will receive the kit bag **FREE**.

**RAPID
INJURY
RELIEF**



Example only



Body Segment	Model No.
Wrist	CRYO/W
Calf	CRYO/C
Ankle	CRYO/A
Elbow	CRYO/E
Knee	CRYO/K
Thigh	CRYO/T
Shoulder	CRYO/S
Back	CRYO/B
Cool Bag	CRYO/Bag
Full Set (8 pieces plus Cool Bag)	CRYO/Set

For more information on other Talar Made products and services visit www.talarmade.com



Order CCT today on
Tel +44 (0)1246 268456
Fax +44 (0)1246 268648
Web www.talarmade.com

TalarMade
 Protect, Support, Correct

CRY010